

CERKAMED

medical company



PRODUCT CATALOG



btr^{pen} Broken
Tool
Remover

ø 0,3 mm



ø 0,5 mm



ø 0,4 mm





CERKAMED

MEDICAL COMPANY



We are **Polish producer** of dental materials



Products sales to over

100 countries!



Cerkamed Team

Quality and **unique design** make our products special



Our medical devices meet the requirements of international standards **ISO 13485:2016**

We sell our products to over 100 countries all over the world, i.a in **Europe, Asia, Africa, and South America**



CE

btrPen




Broken
Tool
Remover

Dental tool for removing broken canal instruments

- + life-long warranty
- + single-tool system for every type of root canal
- + precise grip of the broken instrument
- + thinnest working tip on the market - only \varnothing 0,3 mm



THREE TIPS
ONE GOAL



Set includes:

- BTR PEN – 1 pc stainless steel tool
- BTR CALIBRATOR – 1 pc
- BTR TIP \varnothing 0,3 – 5 pcs
- BTR SPRING – 2 pcs soft and hard (third middle spring fixed on the device)
- Silicone BTR O-RING – 3 pcs

Changeable tips available separately, packed in 5 pcs



Seize the impossible!

L I M I T E D



Golden
edition

Join the world's endodontic elite
with limited edition of BTR PEN

Set includes:

- BTR PEN Limited Golden Edition – 1 pc stainless steel tool
- BTR CALIBRATOR – 1 pc
- GOLDEN BTR TIP – 0,3 mm – 20 pcs
- BTR TIP – 0,4 mm – 10 pcs
- BTR TIP – 0,5 mm – 10 pcs
- BTR SPRING – 2 pcs soft and hard (third middle spring fixed on the device)
- Silicone BTR O-RING – 3 pcs

Let the *luxury* enter your dental office



CE



Ideal kit for training

the removal of broken endodontic files from root canals

1. DIGITAL MICROSCOPE ADSM201

- MINI HDMI input
- MINI USB port
- MICRO SD card input
- hight regulation
- LED lightning

2. MICROSCOPE'S INSTRUCTION FOR USE

3. TWO TRAINING TEETH with pieces of broken endodontic files

HARD 3 deeply placed, blocked pieces of files (requiring processing before removal)

EASY 3 loosely placed pieces of files



- 4. SILICONE TRAINING MASS**
 - transparent and backlit
 - 8 pieces of broken endodontic files
 - a space to place the training tooth
- 5. SILICONE MASS STAND**
 - allows 3D rotation
- 6. BTR Pen WORK METHODOLOGY**
 - work with BTR Pen step by step
- 7. ANTI-SLIP MAT**
 - stability during the work with the kit

BTR DENTAL TRAINING KIT

Also available:



BTR DENTAL TRAINING GUM

- 1. SILICONE TRAINING MASS** with pieces of broken endodontic files
- 2. SILICONE TRAINING MASS STAND**
 - allows 3D rotation
- 3. TWO TRAINING TEETH** with pieces of broken endodontic files (easy/hard)
- 4. ANTI-SLIP MAT**

• CAVIPACK •

Temporary filling

NOVELTY



CAVIPACK

Zinc Oxide, Zinc Sulfate, Calcium Sulfate, excipients.

CE 2409

Temporary filling
3 x 10g

Available soon!

- ready-to-use, zinc oxide and zinc sulfate paste for temporary filling
- it does not contain eugenol
- binds chemically when exposed to moisture
- easy to remove
- radiopaque
- pleasant, rose scent
- economy package **3x10g**



RAINBOW FLOW dentina

LIGHT-CURING FLOW COMPOSITE

- ✓ The RAINBOW FLOW dentine composite is intended to:
 - fill the cavities in deciduous and permanent teeth, in class I, II, III, IV and V class cavities according to Black classification
 - sealing fissures.
- ✓ RAINBOW FLOW dentine lateral teeth in their natural colour.

A2

filler content 62%

RAINBOW FLOW CE 2409
Methacrylate resins

Light-curing flow composite

Available packages:
- syringe 1 g (chosen colour)
+ applicators

- ✓ Package contains:
1 x syringe containing 1 g of preparation + applicators

Application:



✓ Pit and fissure sealing



✓ Filling of permanent teeth and milk teeth or protection of root canal orifices

RAINBOW FLOW dentina
compliant with EN ISO 4049



+ Compression strength



+ Bending strength

RAINBOW FLOW CE 2409
Methacrylate resins
SUPER SIX

Light-curing flow composite

Available packages:
- 6 x syringe 1 g
(colour: orange, blue, pink, yellow,
green, purple) + applicators



- orange
- blue
- pink
- yellow
- green
- purple

3D Stickers



Colour palette for kids

RAINBOW FLOW CE 2409
Methacrylate resins

Light-curing flow composite

Available packages:
- syringe 1 g (chosen colour)
+ applicators



1g

CHLORAXID

Sodium hypochlorite 2% or 5,25% in liquid

CE 2409

PREPARATION for root canals rinsing

Available packages:
CHLORAXID 2% 200 g
CHLORAXID 2% 400 g
CHLORAXID 5,25% 200 g
CHLORAXID 5,25% 400 g



! use with: ENDO-PREP CREAM

ENDO-SOLUTION

EDTA 17% in liquid

CE 2409

LIQUID for root canals widening

Available package:
ENDO-SOLUTION 50 ml



CHLORAXID EXTRA

Sodium hypochlorite 2% or 5,25% in liquid + surfactants

CE 2409

PREPARATION for root canals rinsing

Available packages:
CHLORAXID 2% EXTRA 200 g
CHLORAXID 2% EXTRA 400 g
CHLORAXID 5,25% EXTRA 200 g
CHLORAXID 5,25% EXTRA 400 g



! use with: ENDO-PREP CREAM

CITRIC ACID 40%

40% citric acid solution

CE 2409

LIQUID for root canals rinsing

Available packages:
CITRIC ACID 40% 200 g
CITRIC ACID 40% 400 g



ENDO-PREP CREAM

EDTA 15% in cream, urea peroxide

CE 2409

CREAM for root canals preparation

Available packages:
ENDO-PREP CREAM 2 ml
ENDO-PREP CREAM 5 ml
ENDO-PREP CREAM 10 ml
ENDO-PREP CREAM MEGA PACK 3 x 10 ml



! recommended to use with liquids from CHLORAXID line

ISOPROPYL ALCOHOL

Available packages:
ISOPROPYL ALCOHOL 200 g



ENDO-PREP GEL

EDTA 17% in gel

CE 2409

GEL for root canals preparation

Available packages:
ENDO-PREP GEL 2 ml
ENDO-PREP GEL 5 ml
ENDO-PREP GEL 10 ml





ENDO-TOP **CE**
2409 Endo irrigation needles

Available packages:
ENDO-TOP 100 pcs



ENDO-PACK LABELLED luer - lock dispensers of 5 ml - 20 pcs

Available for products:
CHLORAXID 2% EXTRA,
CHLORAXID 5,25% EXTRA,
CHLORAXID 2%,
CHLORAXID 5,25%,
ENDO-SOLUTION,
CITRIC ACID 40%




20 x



CANAL DETECTOR **CE**
2409 Root canals orifices indicator

methylene blue

Available packages:
CANAL DETECTOR 2 ml




Calcium hydroxide paste for temporary root canal filling

Available packages:
CALCIPAST 2,1 g
CALCIPAST MEGA PACK 4 x 2,1 g



CALCIPAST **CE**
Calcium hydroxide 2409



CALCIUM HYDROXIDE
Calcium hydroxide pure p.a.

Available packages:
CALCIUM HYDROXIDE 10 g



Available packages:
IODOFORM 30 g



IODOFORM



ZINC OXIDE FAST **CE**
2409 Zinc oxide, excipients


Fast-setting root canal repair material

Available packages:
ZINC OXIDE FAST 50 g




Material for making pastes

Available packages:
ZINC OXIDE 50 g



ZINC OXIDE **CE**
Zinc oxide 2409



Oil for making pastes

Available packages:
EUGENOL 10 ml
EUGENOL 20 ml



EUGENOL **CE**
Eugenol oil 100% 2409

• short setting time • plasticity • biocompatibility • strength

BIO MTA+

CE 2409

Biomaterial for filling and rebuilding of root canals

Available packages:
MINI 3 x 0,14 g BIO MTA+ powder, 1 ml liquid
MAXI 10 x 0,14 g BIO MTA+ powder, 2 x 1 ml liquid



MTA+

CE 2409

Material for filling and rebuilding of root canals

Available packages:
MINI 3 x 0,14 g MTA+ powder, 1 ml liquid
MAXI 10 x 0,14 g MTA+ powder, 2 x 1 ml liquid



MTA+ APPLICATOR

CE

Reusable dental tool for handling and application of prepared measure of MTA+ material

Possible variants:
APPLICATOR FOR MTA+ 0,8 mm
APPLICATOR FOR MTA+ 1,2 mm



TOOLS FOR HANDLING AND APPLICATION OF MTA+

CE

Reusable dental tools for handling and application of prepared measure of MTA+ material



MATRIX MTA+

CE

Dosing Block MTA+



MTA SMART SYSTEM
Training set

Endodontic training set for application MTA

Set includes:

MTA SMART SYSTEM TRAINING SET



SILICONE TRAINING MASS

TRAINING TEETH 5 PCS

SILICONE TRAINING MASS STAND
ALLOWS 3D ROTATION

TRAINING MTA POWDER 5 g

TRAINING MTA LIQUID 3 ml

MEASURING CUP FOR MTA POWDER

PIPETTE FOR MTA LIQUID APPLICATION

SLIDE



TRAINING TEETH

5 variants available

<p>2633-101-PCL Maxillary first molar 3 canals, 3 roots Perforation of the bottom of the chamber chamber floor, lateral wall perforation of the palatal root</p>	<p>2111-1A2-PCL Maxillary central incisor 1 canal, 1 root Incomplete root development</p>	<p>2633-102-PCL Maxillary first molar, 3 canals, 3 roots, perforation of the pulp chamber floor, external apex resorption</p>
<p>3623-101-PCL Mandibular first molar 3 canals, 2 roots perforation of the pulp chamber floor</p>	<p>3111-1A1-PCL Mandibular central incisor 1 canal, 1 root, anterior root wall perforation</p>	



Also sell separately!

SYNTEX

CE 2409 Root canals filling and sealing material

Available packages:
syringe of 10 g of the preparation
+mixing tip
+mixing papers



GUTTA PERCHA POINTS

CE 2409 Hand rolled gutta percha points

Available packages:
GUTTA PERCHA 2% 120 pcs
GUTTA PERCHA 2% 120 pcs marked
GUTTA PERCHA 4%, 6% 60 pcs
GUTTA PERCHA 4%, 6% 60 pcs marked
GUTTA PERCHA Pro Taper 60 pcs
GUTTA PERCHA Pro Taper 60 pcs marked
GUTTA PERCHA 2% ROUND PACKAGE 240 pcs
GUTTA PERCHA 2% ROUND PACKAGE 240 pcs. marked
GUTTA PERCHA 4%, 6% ROUND PACKAGE 240 pcs
GUTTA PERCHA 4%, 6% ROUND PACKAGE 240 pcs marked

Color chart and sizes:

08	10	15	20	25	30	35	40	45	50	55	60	65	70	75	80	85	90	95	100
10	15	20	25	30	35	40	45	50	55	60	65	70	75	80	85	90	95	100	105



GUTTA PERCHA PLUGGER

CE Reusable dental tool for vertical condensation of gutta-percha

Available packages:
GUTTA PERCHA PLUGGER 1 pc
GUTTA PERCHA PLUGGER 4 pcs



Oil for dissolving gutta-percha

Available packages:
EUCALYPTOL 10 ml
pipettes for application 5 pcs

EUCALYPTOL

Eucalyptus oil 100%

CE 2409



Oil for dissolving gutta-percha

Available packages:
ORANGE GUTTANE 10 ml
pipettes for application 5 pcs

ORANGE GUTTANE

Orange oil 100%

CE 2409



Available packages:
METHYLENE CHLORIDE 10 ml

METHYLENE CHLORIDE

Dichloromethane pure p.a.



Hand rolled absorbent paper points

Available packages:
ABSORBENT PAPER POINTS 2% 200 pcs
ABSORBENT PAPER POINTS 2% 200 pcs marked
ABSORBENT PAPER POINTS 4%, 6% 100 pcs
ABSORBENT PAPER POINTS 4%, 6% 100 pcs marked
ABSORBENT PAPER POINTS ProTaper 100 pcs
ABSORBENT PAPER POINTS ProTaper 100 pcs marked
ABSORBENT PAPER POINTS 2%
Round package 400 pcs
ABSORBENT PAPER POINTS 2%
Round package 400 pcs marked
ABSORBENT PAPER POINTS 4%, 6%
Round package 240 pcs
ABSORBENT PAPER POINTS 4%, 6%
Round package 240 pcs marked 240 pcs

ABSORBENT PAPER POINTS

CE 2409



ENDO-ASPIRATOR

for removal liquids from the root canal by aspirating

ENDO aspirator^{PRO}

for removal liquids from the root canals

Possible variants:
 ENDO-ASPIRATOR PRO
 ENDO-ASPIRATOR PRO Golden Version
PACKAGE INCLUDES:
 HANDLE - 1 pc
 ELASTIC APPLICATOR - 10 pcs
 CONNECTOR - 3 pcs
 LUER LOCK HOSE - 5 pcs
 MAGNET - 5 pcs
 RAMROD FOR CLEANING THE ENDO-ASPIRATOR PRO - 1 pc



ADDITIONAL ACCESSORIES

PACKAGE INCLUDES:
 LUER LOCK HOSE - 5 pcs
 CONNECTOR - 1 pc




Dental etchant

Available packages:
 BLUE ETCH 2 ml
 BLUE ETCH 10 ml
 BLUE ETCH 50 ml
 BLUE ETCH MEGA PACK 3 x 10 ml
 BLUE ETCH MONSTER PACK 10 x 2 ml




Flowable dental etchant

Available packages:
 BLUE ETCH FLOW 2 ml
 BLUE ETCH FLOW 10 ml
 BLUE ETCH FLOW 50 ml
 BLUE ETCH FLOW MONSTER PACK 10 x 2 ml



Available packages:
 Syringe 2 ml



Available packages:
 Syringe 2 ml



ALUSTAT FOAM CE 2409
Aluminium chloride 20%



FOAM for staunching the bleeding
Available packages:
ALUSTAT FOAM 0,8 g
ALUSTAT FOAM MEGA PACK 4 x 0,8 g



ALUSTAT CE 2409
Aluminium chloride 25%



LIQUID for staunching the bleeding
Available packages:
ALUSTAT 10 g

ALUSTAT GEL CE 2409
Aluminium chloride 25%

GEL for staunching the bleeding
Available packages:
ALUSTAT GEL 5 ml
ALUSTAT GEL 10 ml
ALUSTAT GEL MEGA PACK 3 x 10 ml



TCA DENT CE
Trichloroacetic acid 35%

Preparation for local cauterisation the mucous membranes
Available packages:
dropping bottle 10 ml



BEST-CORD NANO CE 2409
Polyester/Polyamide yarn conjugated. Aluminium chloride impregnated: 000.0/000/00/0: 0,5 ± 0,3 mg/inch 1/2: 1,0 ± 0,4 mg/inch

Knitted retraction cords, 254 cm length




unique smooth weave

- 6 SIZES:
- 000.0
 - 000
 - 00
 - 0
 - 1
 - 2



RETRACTION CORD PACKERS CE



Tools for packing retraction cords in the gingival pockets
Possible variants:
RETRACTION CORD PACKERS
1. UNIVERSAL
2. ROUND 1,5 mm
3. ROUND 3 mm



Each item of Rubber dam Set available

separately!


12 types of clamps
Also available in satin version

RUBBER-DAM SET



Tools for rubber dam

- Set includes:
RUBBER-DAM PUNCH
CLAMPS FORCEPS
RUBBER-DAM FRAMES (3 sizes)
THE SET OF CLAMPS
BASE FOR CLAMPS



Flat jaws for small lower molars (No 3).
Festooned wings for large upper molars.
Wings fit the gum curvature (No 5).
Festooned jaws for molar's roots (No W8A).
Serrated jaws for lower right and upper left molars (No12A).
Serrated jaws for lower left and upper right molars (No 13A).
Small and large jaw for anteriors (No 210).
Flat jaws for anteriors (No 212).
Rounded jaws for premolars (No 27).
Rounded and at jaws for premolars (No 29).
Small jaws and high bow for anteriors and premolars (No 00).
Wings, flat and large jaws for premolars (No 2A).
Festooned wings, flat jaws for large premolars (No 208).

ADDITIONALLY RECEIVE for free!

- RUBBER-DAM (medium plain)
- RUBBER-DAM TEMPLATE PREMIUM with pen
- RUBBER-DAM LIQUID
- Light-cured rubber dam liquid
- RUBBER-DAM NAPKINS
- RUBBER-DAM CORD
- Dental dam stabilizing cord (medium)

RUBBER-DAM SET MINI



Tools for rubber dam Start Pack

- Set includes:
RUBBER DAM PUNCH
CLAMPS FORCEPS
RUBBER DAM FRAME (M SIZE)
CLAMPS NO 12A, 13A, 208, 210



Light-Cured RUBBER-DAM liquid

Available packages:
RUBBER-DAM LIQUID 1,2 ml
RUBBER-DAM LIQUID MEGA PACK 4 x 1,2 ml

RUBBER-DAM LIQUID



urethane resins,
inorganic fillers



Available package:
36 pcs
Available variants:
- plain BLUE: medium 0,18 mm
or heavy 0,22 mm
- mint GREEN: medium 0,18 mm
or heavy 0,22 mm
- mint BLACK: medium 0,18 mm

RUBBER-DAM



Available packages:
RUBBER-DAM NAPKINS 36 pcs

RUBBER-DAM NAPKINS



Dental dam stabilizing cord

Available packages:
RUBBER-DAM CORD fine
RUBBER-DAM CORD medium
RUBBER-DAM CORD thick

RUBBER-DAM CORD

Package containing 1,2 m of rubber dam cord

Available colours:

- black
- purple
- blue
- red
- green
- yellow
- white



Silicone rings for marking of instruments in stomatology

Available packages:
SILICONE RINGS 20 pcs
SILICONE RINGS 100 pcs

SILICONE RINGS



GLUCOSITE

aqua, chlorhexidine digluconate, alcohol, hydrogen peroxide

LIQUID for gingival pockets

Available packages:
GLUCOSITE 2 ml
GLUCOSITE 50 ml
GLUCOSITE MONSTER PACK
10 x 2 ml

GEL for gingival pockets

Available packages:
GLUCOSITE GEL 2 ml
GLUCOSITE GEL 50 ml
GLUCOSITE GEL MONSTER PACK
10 x 2 ml

GLUCOSITE GEL

aqua, glycerine, chlorhexidine digluconate, methylcellulose, alcohol, hydrogen peroxide

INOX

propanotriol (glycerol), Ph Eur purified water, gel base



STOP
Inhibition

INOX

Gel promoting complete polymerisation of dental composites

Available packages:
INOX 2 ml
INPX MEGA PACK 4 x 2 ml

RED DETECTOR

Rhodamine



2409

Caries detector

Available packages:
RED DETECTOR 2 ml

PULP SPRAY



SPRAY for pulp vitality testing

Available packages:
PULP SPRAY 200 ml

Aroma:
- mint
- orange

elastic applicators


for application



Elastic applicators for liquids, creams, pastes, gels

Available packages:
20 pcs

APPLICATORS




◀ ELASTIC APPLICATOR 27G straight
length: 41 mm
diameter: 0,4 - 0,2 mm




◀ ELASTIC APPLICATOR 27G bent
length: 41 mm
diameter: 0,4 - 0,2 mm




◀ ELASTIC APPLICATOR 25G bent
length: 41 mm
diameter: 0,5 - 0,3 mm



◀ ELASTIC APPLICATOR 25G long bent
length: 56 mm
diameter: 0,5 - 0,3 mm



◀ ELASTIC APPLICATOR 22G curve
length: 32,3 mm
diameter: 0,7 - 0,5 mm



◀ ELASTIC APPLICATOR 21G long bent
length: 56 mm
diameter: 0,8 - 0,6 mm



◀ ELASTIC APPLICATOR 27G bent with two side holes
length: 35 mm
diameter: 0,4 mm

ADDITIONAL ACCESSORIES

Available packages:
1. ELASTIC APPLICATORS
10 pcs
ELASTIC APPLICATORS
20 pcs

2. GEL APPLICATORS
20 pcs
GEL APPLICATORS
50 pcs

3. BRUSH APPLICATORS
10 pcs
BRUSH APPLICATORS
20 pcs

4. APPLICATORS FOR ALLUSTAT FOAM
10 pcs

5. APPLICATORS FOR SYNTAX
5 pcs

6. PAPER FOR PRODUCT MIXING
30 pcs

7. APPLICATION PIPETTES
5 pcs



for suction



PERFECT PROTECTION

PERFECT PROTECTION

CE

Available colours:

- green
- blue
- pink
- beige
- white
- transparent

COMFORT
with set of foils



COMFORT LIGHT
with set of foils



SET of removable protective foils

Available packages:
4 x RELAX (0,175 mm)
5 x STANDARD (0,125 mm)
4 x LIGHT RELAX (0,175 mm)

ANTISTATIC FOILS

CE

COMFORT OF WORK

PROTECTIVE VISORS ATTACHED TO GLASSES

CE



ONE SIZE


PROTECTIVE SHIELDS FOR A PATIENTS

CE

Packing includes:
Shield frame with plastic hooks
+ set of exchangeable foils - 8 pcs

Available colours:

- green
- blue
- pink
- beige
- white
- transparent



PROTECTIVE PADS FOR A PATIENT

CE




Packing includes:
20 pcs OF DISPOSABLE PADS



GLASSES FOR POLYMERIZATION

CE





ENDO-STATION CE

Line of endodontic organizers

- Available in colors*:
- green
 - blue
 - transparent

* ENDO-STATION GO GUTTA BOX available only transparent

ENDO-STATION GO MINI CE

- Available in colors:
- green
 - blue
 - transparent

DENTAL CONTRASTOR CE

- universal
- single use
- flexible



Protective enamel with fluoride for children

Available packs:
 - FLUOR DEFENDER MINI
 3 x 1 ml of preparation
 + set of applicators
 - FLUOR DEFENDER MAXI
 10 x 1 ml of preparation
 + set of applicators

FLUOR DEFENDER CE

0,1% fluorosilan (1600 ppm F)
 HEMA

MARKETING MATERIALS FOR PATIENT



- Available packages:
1. 3D STICKERS 10 pcs, 50 pcs
 2. LEAFLET „ROOT CANAL TREATMENT“ DEDICATED FOR PATIENT 20 pcs, 50 pcs, 100 pcs
 3. PUPPET 1 pcs, 2 pcs, 3 pcs
 4. BOXES FOR MILK TEETH 5 pcs
 5. DIPLOMA „BRAVE PATIENT“ mix 20 pcs
 6. GUESSING GAMES FOR A BRAVE PATIENT 20 pcs
 7. TEETH BRUSHING CALENDAR for kids
 Set included: teeth brushing calendar, diploma, stickers, magnet
 8. PUZZLES FOR CHILDREN 36 ELEMENTS

Perfect picture!



15
 pcs
 Anti-reflex



Specialised Mouthwash

- Available packages:
- Start Pack concentrate in syringe 50 ml + bottle + dispenser
 - concentrate in syringe 50 ml

- Two variants available:
- MINT
 - ORIGINAL



6in1



JOIN TO CERKAMED IN:



PPH CERKAMED

Wojciech Pawłowski



Address: ul. Kwiatkowskiego 1
37-450 Stalowa Wola
Poland



Phone: +48 15 840 08 09
+48 15 842 35 85, int. 809
+48 785 665 505



E-mail: export@cerkamed.pl
WWW: www.cerkamed.com